

BioDrive **Cannulated** Screws Titanium or Stainless Steel?



The Choice Is **Yours**

BOMET[®]

The versatile Biodrive Cannulated Screw system has been developed to address a multi-application osteosynthesis. The range of sterile and non-sterile Stainless Steel screws in diameter 2.8, 3.5, 4.0, 4.5, 5.0, 6.5 and 7.5mm is now completed by Sterile partially threaded 4.0mm and 6.5mm diameter Titanium screws.

Cannulated design

For an easy insertion and accurate placement.
Screws can be implanted percutaneously

Low profile head easily countersunk

Reduces bone damage and soft tissue irritation

Self-Drilling & Self-Tapping

Simplifies the procedure, shorten the operative time

Large diameter guide wires

Reduce the chance of guide deviation upon insertion



Reverse cutting threads

Facilitates screw removal and may help to reduce screw breakage upon extraction

Precise, color coded instrumentation

Provides clear and efficient use

Titanium Alloy Sterile Screws

High fatigue strength alloy, non-magnetic MRI compatible. Addresses nickel sensitive patients

Ordering Information

Cat. No.	Description	Cat. No.	Description
3212-040-010S	BIODRIVE Ti Ø4.0X10MM Screw (Part Thd)	3213-065-035S	BIODRIVE Ti Ø6.5X35MM Screw (16MM Thd)
3212-040-012S	BIODRIVE Ti Ø4.0X12MM Screw (Part Thd)	3213-065-040S	BIODRIVE Ti Ø6.5X40MM Screw (16MM Thd)
3212-040-014S	BIODRIVE Ti Ø4.0X14MM Screw (Part Thd)	3213-065-045S	BIODRIVE Ti Ø6.5X45MM Screw (16MM Thd)
3212-040-016S	BIODRIVE Ti Ø4.0X16MM Screw (Part Thd)	3213-065-050S	BIODRIVE Ti Ø6.5X50MM Screw (16MM Thd)
3212-040-018S	BIODRIVE Ti Ø4.0X18MM Screw (Part Thd)	3213-065-055S	BIODRIVE Ti Ø6.5X55MM Screw (16MM Thd)
3212-040-020S	BIODRIVE Ti Ø4.0X20MM Screw (Part Thd)	3213-065-060S	BIODRIVE Ti Ø6.5X60MM Screw (16MM Thd)
3212-040-022S	BIODRIVE Ti Ø4.0X22MM Screw (Part Thd)	3213-065-065S	BIODRIVE Ti Ø6.5X65MM Screw (16MM Thd)
3212-040-024S	BIODRIVE Ti Ø4.0X24MM Screw (Part Thd)	3213-065-070S	BIODRIVE Ti Ø6.5X70MM Screw (16MM Thd)
3212-040-026S	BIODRIVE Ti Ø4.0X26MM Screw (Part Thd)	3213-065-075S	BIODRIVE Ti Ø6.5X75MM Screw (16MM Thd)
3212-040-028S	BIODRIVE Ti Ø4.0X28MM Screw (Part Thd)	3213-065-080S	BIODRIVE Ti Ø6.5X80MM Screw (16MM Thd)
3212-040-030S	BIODRIVE Ti Ø4.0X30MM Screw (Part Thd)	3213-065-085S	BIODRIVE Ti Ø6.5X85MM Screw (16MM Thd)
3212-040-032S	BIODRIVE Ti Ø4.0X32MM Screw (Part Thd)	3213-065-090S	BIODRIVE Ti Ø6.5X90MM Screw (16MM Thd)
3212-040-034S	BIODRIVE Ti Ø4.0X34MM Screw (Part Thd)	3213-065-095S	BIODRIVE Ti Ø6.5X95MM Screw (16MM Thd)
3212-040-036S	BIODRIVE Ti Ø4.0X36MM Screw (Part Thd)	3213-065-100S	BIODRIVE Ti Ø6.5X100MM Screw (16MM Thd)
3212-040-038S	BIODRIVE Ti Ø4.0X38MM Screw (Part Thd)	3213-065-105S	BIODRIVE Ti Ø6.5X105MM Screw (16MM Thd)
3212-040-040S	BIODRIVE Ti Ø4.0X40MM Screw (Part Thd)	3213-065-110S	BIODRIVE Ti Ø6.5X110MM Screw (16MM Thd)
3212-040-042S	BIODRIVE Ti Ø4.0X42MM Screw (Part Thd)	3213-065-115S	BIODRIVE Ti Ø6.5X115MM Screw (16MM Thd)
3212-040-044S	BIODRIVE Ti Ø4.0X44MM Screw (Part Thd)	3213-065-120S	BIODRIVE Ti Ø6.5X120MM Screw (16MM Thd)
3212-040-046S	BIODRIVE Ti Ø4.0X46MM Screw (Part Thd)	3213-065-125S	BIODRIVE Ti Ø6.5X125MM Screw (16MM Thd)
3212-040-048S	BIODRIVE Ti Ø4.0X48MM Screw (Part Thd)	3213-065-130S	BIODRIVE Ti Ø6.5X130MM Screw (16MM Thd)
3212-040-050S	BIODRIVE Ti Ø4.0X50MM Screw (Part Thd)	191151S	BIODRIVE WASHER TITANIUM SMALL
3213-065-030S	BIODRIVE Ti Ø6.5X30MM Screw (16MM Thd)	191153S	BIODRIVE WASHER TITANIUM LARGE